

## Simultaneous detection and genotyping of 35 HPV genotypes



- ✓ Liquid- based cytology
- ✓ Cytological and anal swabs
- ✓ Paraffin- embedded samples



No DNA  
extraction

### High Risk:

16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56,  
58, 59, 66, 68, 73 and 82

### Low Risk:

6, 11, 40, 42, 43, 44, 54, 55, 61, 62, 67, 69, 70,  
71, 72, 81 and 84

<b>B</b>	33	58	42	71	16	52	<b>B</b>	
<b>B</b>	35	59	43	72	18	53	6	69
<b>C</b>	39	66	44/55		26	56	11	70
<b>U</b>	45	68	54	84	31	58	40	71
16	51	73	61	<b>B</b>	33	59	44/55	72
18	52	82	62/81	<b>C</b>	35	66	54	
26	53	6	67	<b>U</b>	39	68	61	84
31	56	11	69	42	45	73	62/81	
	<b>B</b>	40	70	43	51	82	67	