

Simultaneous detection of 12 STD-related pathogens



- ✓ Urine
- ✓ Semen
- ✓ Endocervical swabs
- ✓ Pharyngeal swabs
- ✓ Perianal swabs
- ✓ Urethral swabs
- ✓ Liquid cytology
- ✓ Vaginal swabs



No DNA
extraction

Organism

- _____ Mycoplasma genitalium
- _____ Mycoplasma hominis
- _____ Trichomonas vaginalis
- _____ Ureaplasmas (U. urealyticum / U. parvum)
- _____ Haemophilus ducreyi
- _____ Herpesvirus simplex type 1/ Herpesvirus simplex type 2
- _____ Herpesvirus simplex type 2
- _____ Treponema pallidum
- _____ Neisseria gonorrhoeae
- _____ Chlamydia trachomatis (Serovares A-K)
- _____ Chlamydia trachomatis (Serovares A-L3)

B	MG				CTS-1		B	
B		MH				CTS-2		
CI	TV		UU-P		MG		UU-P	
BG		HD				MH		
	HSV-1/ HSV-2		NG	B	TV		NG	
		HSV-1		CI		HD		
	CTS-1		TP	BG	HSV-1/ HSV-2		TP	
		CTS-2				HSV-1		
	B							